6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Certificate of Analysis

Product: Purified Human IgA1

Part Number: Hu-107-M Lot Number: 80-101-052621

Concentration: $2.81 \text{ mg/ml} (E^{1\%} \text{ at } 280 \text{ nm} = 13.0)$

Amount: 1.0 mg

Form: Clear, colorless liquid, 0.2 µm filtered

Purity: > 95 % based on SDS-PAGE

Isotype: Kappa

Buffer: 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2

Preservative: 0.05 % (w/v) Sodium Azide

Storage: -20° C for long term storage. 2-8° C for short term storage.

Expiration Date: May 26, 2024

Reactivity: Single arc by IEP against antisera to whole serum, human IgA1, and human

kappa. No reaction to antisera to human IgA2, IgG, IgD, IgE, and lambda.

Country of Origin: Human myeloma plasma was obtained from donors of US origin.

Disease Testing: Donor serum was tested by FDA approved methods and found negative for

antibodies to HIV 1/2 and HCV. Donor serum was tested by FDA approved

methods and found negative for HBsAg, HIV-1.

Applications: This antibody is suitable for use in most immunoassay formats as a control or

standard. The optimal working dilution should be determined by the

investigator.

Disclaimer: For in vitro Laboratory Use Only. Not for diagnostic or therapeutic use. Not

for human or animal consumption. The proper selection and use of our products is the sole responsibility of the end user, and therefore we can offer no guarantee to a specific experimental outcome. Suggested applications of our products are not recommendations to use our products in violation of any

patent or as a license under any patent of ImmunoReagents, Inc.

Product may not be resold or modified for resale without prior written

approval of ImmunoReagents, Inc.

Phone: 919-831-2240 · www.ImmunoReagents.com · Fax: 919-831-2240

6003 Chapel Hill Road · Suite # 153 · Raleigh, NC 27607

Testing Results

PRODUCT: Purified Human IgA1

PART NUMBER: Hu-107-M **LOT NUMBER:** 80-101-052621

M

<u>TESTS</u> <u>SPECIFICATIONS</u> <u>RESULTS</u>

Appearance: Clear, colorless liquid PASS

Purity: Two predominant bands at 55 and 24 kDa based PASS


on SDS-PAGE

Reactivity: Single arc by IEP against antisera to whole serum, PASS

human IgA1, and human kappa. No reaction to antisera to human IgA2, IgG, IgD, IgE, and lambda.

SDS-PAGE

1 2



- 1) Hu-107-M, Lot # 80-101-052621
- 2) Mol Weight Standard

Date: May 27, 2021 Signature: Quality Control

Phone: 919-831-2240 · www.ImmunoReagents.com · Fax: 919-831-2240



Phone: 919-831-2240 www.ImmunoReagents.com Fax: 919-831-2240